Size = A3 (420X297mm), Black Print

# TrueSafe® Safety I.V. Catheter English

NOTE: U.S. Federal Law restricts this device sale by or on order of registered physician.


INTENDED USE: The MedSource TrueSafe Safety I.V. Catheter is indicated to sample blood or administer fluids intravenously. The MedSource TrueSafe® Safety I.V. Catheter may be used for any patient population with consideration given to adequacy for the vascular anatomy appropriateness for the solution being administered and duration of the therapy.

The following instructions contain important information regarding the safe use of this product. Please read the Instruction For Use contents entirely before using this product. Failure to properly follow warnings and instructions could result in serious injury or death to the patient and/or clinician.

### Product Description:

The device TrueSafe® Safety I.V. Catheter consists of major components as in fig. (1), fig. (2) & fig. (3), (A) Stainless Steel Needle, (B) Front chamber, (C) Gauge Chamber, (D) Catheter Cannula, (E) Safety Cap, (F) Stainless Steel Spring, (G) Flash Chamber, (H) Activation Button, (I) Catheter Hub.

The TrueSafe® Safety I.V. Catheter gauges are color coded for positive identification and range from 14G to 24G  
ORANGE = 14G GREY = 16G GREEN = 18G PINK = 20G BLUE = 22G YELLOW = 24G

### Warnings:

Clinicians must follow Instructions for Use. In addition to following the Instructions for use, clinicians must be trained to ensure proper usage of the device. Failure to do either of these could result in death or serious injury to the patient and/or clinician. To avoid an accidental needle stick, use the device as described in the Instructions for Use and discard immediately into an approved sharps container. Do not apply pressure during usage that causes the needle / catheter to flex. In the unlikely event that the device does not function as described in the Instructions for Use, immediately dispose of the device in an approved sharps container. Prior to insertion, caution needs to be taken to ensure that the catheter does not come in contact with any foreign objects. Contact with a foreign object creates the potential for embolus. Never reinsert the introducer needle inside the catheter. Reinsertion threatens the integrity of the catheter and could lead to embolus.

### Cautions:

This device is in an unopened package and is sterile, non-toxic and non-pyrogenic. Discard if open, wet or damaged. The materials used to manufacture this device Not made with natural rubber latex. This product is not made with PVC. Not made with DEHP. This device is Radiopaque.

This device is for single use only. Any attempt to re-sterilize and/or recondition the device will compromise the integrity of the device and may have an adverse affect on performance. US Federal Law restricts the sale of this device to or on the order of a registered physician.

### Precautions:

This device is designed to reduce needle stick occurrences. Care must always be exercised when using this product. In addition to the Instructions for Use provided, it is strongly recommended that the clinician follow the guidelines put forth by the CDC and OSHA standards (USA) for bloodborne pathogens when using this product.

During catheter insertion, maintenance and removal follow current "Guidelines for the Prevention of Intravascular Catheter-Related Infections," Centers for Disease Control and Prevention, Atlanta, GA (USA) Institutional policy and procedures, and professional standard of care.

Securely affix the catheter to the patient. Improper securing can result in loss of vascular access. A firm connection between the catheter hub and administration set is necessary to prevent leaking. The introducer needle is longer than the catheter tube and will enter the vein or artery first. This can cause a blood flashback prior to the entrance of the catheter tube into the vein or artery. Advancing the catheter and needle together slightly will ensure intravascular access for the catheter. When advancing slightly, it is recommended to reduce the angle of approach in order to minimize the chance of puncturing the posterior wall of the vessel.

In the event of an unsuccessful venipuncture, stabilize the catheter hub and activate the safety mechanism. The safety position can be found by hearing a slight "click" and noticing the retraction of the needle into the gauge chamber via the spring. Remove the catheter and properly dispose of both the spent catheter and locked Safety I.V. Catheter assembly.

### Instruction for use:

Due to the risk of bloodborne pathogen exposure, follow Universal Precautions during placement, use and removal of TrueSafe Safety I.V. Catheter.

Clinicians must be trained in the practice of venipuncture and be aware of inherent dangers. Keep unsharpened needles away from body and fingers at all times.

Aspectic technique, proper skin preparation and continued protection of the site are essential. Observe universal precautions on all patients.

A. Remove safety cap in a straight outward motion and inspect catheter unit.

B. Grasp device front chamber between thumb and fingers, verify the advancement recto and needle bevel are in the up position.

C. Perform venipuncture.


D. Observe blood return. **Sizes 20, 22, 24, GA: blood return will be visible through the catheter flash window** Fig. (3). Lower and advance the catheter 1/8".

E. Keeping the device stable, advance the catheter off of the hub and into the vein.

F. Before withdrawing the needle from the catheter, depress the activation button to activate the safety mechanism. The needle should retract fully into the gauge chamber.

G. Connect I.V. Administration line.

H. Immediately discard unit into an approved sharps container.

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
Manufactured for MedSource, 8600 Shelby Court Chanhassen, MN 55317, USA  
www.medsourceclabs.com

	Product Ref. No.		Keep dry
	Batch Number		Fragile, handle with care
	Date of Manufacturing		This way up
	Use By		Quantity
	Do Not Reuse		Non-pyrogenic
	Sterilized by Ethylene Oxide Gas		Manufacturer
	Caution, Consult Accompanying Documents		Authorized EC Representative in the European Community
	Consult Instructions for Use		Recyclable Packaging
	Do Not Use if Package is Damaged		Green Dot
	Do not resterilize		Recyclable
	Keep away from sunlight		

**MEDSOURCE**

Manufactured for:  
MedSource  
8600 Shelby Court  
Chanhassen, MN 55317, USA

Rx only

Made in India  
AW/IFU\_BP, Rev. 01  
Issue Date: 24.06.21  
N. Code: HR/DEVICE/MFG/MD/2019/000135

# Size = A3(420X297mm), Black Print

## TrueSafe Comfort™ Safety I.V. Catheter English


**CAUTION:** U.S. Federal Law restricts this device to sale by or on the order of a registered physician.

**INTENDED USE:** The MedSource TrueSafe™ Safety I.V. Catheter is indicated to sample blood or administer fluids intravenously. The MedSource TrueSafe™ Safety I.V. Catheter may be used for any patient population with consideration given to adequacy for the vascular anatomy appropriate for the solution being administered and duration of the therapy.

**! The following instructions contain important information regarding the safe use of this product. Please read the Instruction For Use contents entirely before using this product. Failure to properly follow warnings and instructions could result in serious injury or death to the patient and/or clinician.**

**Product Description:**  
The device TrueSafe™ Safety I.V. Catheter consists of major components as in Fig. 1), Fig. 2) & Fig. 3). (A) Stainless Steel Needle, (B) Front Chamber, (C) Gauge Chamber, (D) Catheter Cannula, (E) Safety Cap, (F) Stainless Steel Spring, (G) Flash Chamber, (H) Activation Button, (I) Catheter Hub.

The TrueSafe™ Safety I.V. Catheter gauges are color coded for positive identification and range from 14G to 26G  
ORANGE = 14G GREY = 16G GREEN = 18G PINK = 20G BLUE = 22G YELLOW = 24G PURPLE = 26G

**Warnings:**  
Clinicians must follow Instructions for Use. In addition to following the Instructions for use, clinicians must be trained to ensure proper usage of the device. Failure to do either of these could result in death or serious injury to the patient and/or clinician.  
To avoid an accidental needle stick, use the device as described in the Instructions for Use and discard immediately after an approved sharps container.  
Do not apply pressure during usage that causes the needle/catheter to flex.  
In the unlikely event that the device does not function as described in the Instructions for Use, immediately dispose of the device in an approved sharps container.  
Prior to insertion, caution needs to be taken to ensure that the catheter does not come in contact with any foreign objects. Contact with a foreign object creates the potential for embolus.  
Never insert the introducer needle into the catheter. Reinsertion threatens the integrity of the catheter and could lead to embolus. TrueSafe™ is intended for short-term use, 72 hours. Usage beyond this lead to infection.


**Precautions:**  
This device is in an unopened package and is sterile, non-toxic and non-pyrogenic. Discard if open, wet or damaged.  
The materials used to manufacture this device do not contain natural rubber latex.  
This product is not made with PVC. This product is not made with DEHP. This device is Radioopaque.

This device is for single use only. Any attempt to re-sterilize and/or recondition the device will compromise the integrity of the device and may have an adverse affect on performance.  
**CAUTION:** U.S. Federal Law restricts this device to sale by or on the order of a registered physician.

**Precautions:**  
This device is designed to reduce needle stick occurrences. Care must always be exercised when using this product. In addition to the Instructions for Use provided, it is strongly recommended that the clinician follow the guidelines put forth by the CDC and OSHA standards (USA) for bloodborne pathogens when using this product.  
During catheter insertion, maintenance and removal follow current Guidelines for the Prevention of Intravascular Catheter-Related Infections, Centers for Disease Control and Prevention, Atlanta, GA (USA) Institutional policy and procedures, and professional standard of care.  
Securely affix the catheter to the patient. Improper securing can result in loss of vascular access.  
A firm connection between the catheter hub and administration set is necessary to prevent leaking.  
The introducer needle is longer than the catheter tube and will enter the vein or artery first. This can cause a blood flashback prior to the entrance of the catheter tube into the vein or artery. Advancing the catheter and needle together slightly will ensure intravascular access for the catheter. When advancing slightly, it is recommended to reduce the angle of approach in order to minimize the chance of puncturing the posterior wall of the vessel.  
In the event of an unsuccessful venipuncture, stabilize the catheter hub and activate the safety mechanism. The safety position can be found by hearing a slight "click" and noticing the retraction of the needle into the gauge chamber via the spring. Remove the catheter and properly dispose of both the spent catheter and locked Safety I.V. Catheter assembly.

**Instruction for use:**  
**Due to the risk of bloodborne pathogen exposure, follow Universal Precautions during placement, use and removal of TrueSafe™ Safety I.V. Catheter.**  
Clinicians must be trained in the practice of venipuncture and be aware of inherent dangers.  
Keep unutilized needles away from body and fingers at all times.  
Aseptic technique, proper skin preparation and continued protection of the site are essential.  
Observe universal precautions on all patients.

- Remove safety cap in a straight outward motion and inspect catheter unit.
- Grasp device front chamber between thumb and fingers, verify the advancement tab and needle level are in the up position.
- Perform venipuncture.
- Observe blood return. **Sizes 20, 22, 24, GA: blood return will be visible through the catheter flash window**
- Fig. 3) Lower and advance the catheter 1/8".
- Keeping the device stable, advance the catheter off of the hub and into the vein.
- Before withdrawing the needle from the catheter, depress the activation button to activate the safety mechanism. The needle should retract fully into the gauge chamber.
- Connect I.V. Administration line.
- Immediately discard unit into an approved sharps container.

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
	Product Ref. No.		Fragile, handle with care
	Batch Number		Quantity
	Date Of Manufacturing		Non pyrogenic
	Use By		Manufacturer
	STERILE   EC Sterilized by Ethylene Oxide Gas		Authorized EC Representative in the European Community
	Caution, Consult Accompanying Documents		Recyclable Packaging
	Do Not Use if Package is Damaged		Green Dot
	Do not resterilize		Recycle
	Keep away from sunlight		Rx only
	MedSource Labs 8600 Shelby Court Chanhassen, MN 55317		AWIFU_BQ_Rev. 01 Issue Date: 03.03.19 HR-DRUGS-637-B (H)