

# BAIRD STRAP

## 5 STRAPS PER POUCH

MS-95015

### Directions:

1. Place strap on patient's body so the shoulder straps lie across patient's shoulders and the center strap lies along the patient's anterior midline, hook and loop side up.
2. Starting at one side of the end body strap at the lower leg, thread loop material strap through the handhold on backboard. Pull snugly and press onto hook material strap to secure, repeat on other side.
3. Repeat about process with 4 remaining body straps.
4. Secure the 2 straps on patient's shoulders; thread loop material strap through the handhold on the backboard. Pull snugly and press onto hook material strap to secure. Repeat on other side.
5. Repeat above procedure with 2 straps coming off end strap at patient's feet. Check all straps again.

**To Release:** Begin with shoulder straps and separate hook from loop on each strap. Repeat procedure with all straps.



### WARNING:

Improper use could result in injury or death. Only properly trained healthcare professionals should use the Baird Strap. All persons using the Baird Strap should be thoroughly trained for safe application before using. The instructions are as a guide only. Use your local guidelines, procedures, and medical controls.

LOT



Single Use Only

REV. 02