



[www.medsourcelabs.com](http://www.medsourcelabs.com)

## **Foldable Mask Directions for Use MS-PM103W**

The Foldable Mask provides a physical barrier between the rescuer and victim, eliminating direct contact of the rescuer's lips with the unknown subject.

The Foldable Mask promotes an airtight seal to the face allowing ventilation through both the mouth and nose simultaneously. Proper training in ventilation and CPR should be obtained before using this device. Prior training of the Foldable Mask on manikins is recommended.

### **USE ON ADULTS AND CHILDREN:**

1. Remove the mask and valve assembly from case. Push out the soft dome with your fingers (Fig. 1) and ensure that the one way valve is in place (Fig. 2).
2. Clear the airway of any debris.
3. Place the mask over the patient's mouth and nose (Fig. 3). Use the rim of the mask to retract the lower lip so that the patient's mouth remains open under the mask. Use thumbs to hold the mask and your fingers to grasp the lower jaw in front of the earlobes. Pull upward to open the airway and seal the mask (Fig. 4).
4. Blow slowly into the airway valve twice and watch as the patient's chest rises (Fig. 5).
5. Ventilate the patient every five seconds for an adult or every three seconds for a child. Continue rescue breathing until the patient responds or you are relieved by high medical authority.

**NOTES:** If patient vomits, remove the Foldable Mask and clean both the patient's airway and the mask with your gloved fingers. If you should contaminate either the valve, mouthpiece or any of your exposed skin surfaces, use the alcohol wipes to clean these areas immediately. Resume rescue breathing promptly.


### **USE WITH OXYGEN (OPTIONAL):**


The Foldable Mask may be used with supplemental oxygen. Attach the oxygen tubing to the optional oxygen valve and set the flow rate. A head strap is included on all models that have the optional oxygen valve. If the patient is spontaneously breathing, use the head strap to position the Foldable Mask on the patient's face. If the patient is not breathing, begin mouth to mask ventilation immediately. The use of supplemental oxygen to the Foldable Mask will increase the delivery of oxygen to the patient. The rescuer will not be affected by the use of supplemental oxygen.

**USE WITH OTHER VENTILATION EQUIPMENT:**  
The Foldable Mask may be used with a Bag-Valve Mask or oxygen powered ventilation devices. Simply detach the valve filter assembly and attach the Foldable Mask to the device.








**CLEANING:**  
The Foldable Mask is reusable and easy to clean. Wash the mask in warm soapy water. Rinse in clean water. Use a 1:64 bleach to clean water mixture to soak the mask for ten minutes. Rinse the mask again with fresh water and allow to air dry. The face mask can't be disinfected using an activated glutaraldehyde solution (such as Cidex). The one-way valve and filter housing is intended for single patient use and should be discarded after use.

**DO NOT** pasteurize, boil, steam or gas autoclave this device.





 **WARNING:** This product can expose you to chemicals including Di(2-ethylhexyl)phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).



	Importer		Lot number	 high temp low temp Temperature limit
	Caution		Control/date of manufacture	
	Country of manufacture		Expiration	

**MEDSOURCE**   
[www.medsourcelabs.com](http://www.medsourcelabs.com)

 8600 Shelby Court  
Chanhassen MN, 55317 USA

Not made with natural rubber latex.  
 Made in China  
MSL-REG-LQR-8536, REV. 02