



MedSource
Diagnostics

Sprague Stethoscope

The MedSource Sprague Stethoscope features a double-tube configuration that provides improved acoustics. This set includes an adult and pediatric diaphragm and three bells (1, 1.25, 1.75) for accurate detection of high, middle, and low frequency sounds in adults and children. The stethoscope is 22" long and has an adjustable chrome-plated brass binaural.



Specification:

- Length: 56cm = 22"
- Reversible chestpiece with large and small diaphragm
- (1) Pair of chrome plated earpieces with plastic eartips
- Heavy wall 22" tubing with metal clip
- Interchangeable adult, pediatric and infant bell type chestpieces
- Spare ear tips, large and small sizes
- Spare transparent plastic diaphragms, large and small sizes
- Plastic pouch for accessories
- Double tube configuration for reliable and clear acoustics
- Adjustable high quality, chrome plated brass earpieces
- Comfort sealing soft ear tips
- Durable high quality PVC Tubes

Accessories Included:



Rev. 1

www.gomedsource.com

Item #	Description	UOM
MS-STBK	Sprague Stethoscope, Black	EA, 50 EA/CS
MS-STBLU	Sprague Stethoscope, Blue	EA, 50 EA/CS
MS-STR	Sprague Stethoscope, Red	EA, 50 EA/CS
MS-STG	Sprague Stethoscope, Green	EA, 50 EA/CS
MS-STY	Sprague Stethoscope, Yellow	EA, 50 EA/CS