


ClearSafe®

Slide Style Safety IV Catheter



- PTFE catheter material
- Slide style needle retraction
- Ground stainless-steel needle
- Straight hub
- Disposable
- Non-pyrogenic
- Not made with natural rubber latex
- Not made with DEHP
- ISO compliant color coding
- Needlestick Safety and Prevention Act (NSPA) of 2000 compliant




www.medsourcelabs.com

Item #	Description	UOM
MS-84114	ClearSafe® Safety IV, 14G x 1 1/4"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84116	ClearSafe® Safety IV, 16G x 1 1/4"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84118	ClearSafe® Safety IV, 18G x 1 1/4"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84120	ClearSafe® Safety IV, 20G x 1 1/4"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-841202	ClearSafe® Safety IV, 20G x 1"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84122	ClearSafe® Safety IV, 22G x 1"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84124	ClearSafe® Safety IV, 24G x 3/4"	50 EA/BX, 4 BX/CS, 200 EA/CS

ClearSafe Comfort[®]

Slide Style Safety IV Catheter



SPECIFICATIONS



- PUR catheter material
- Slide style needle retraction
- Ground stainless-steel needle
- Straight hub
- Disposable
- Non-pyrogenic
- Not made with natural rubber latex
- Not made with DEHP
- ISO compliant color coding
- Needlestick Safety and Prevention Act (NSPA) of 2000 compliant



www.medsourcelabs.com

Item #	Description	UOM
MS-84214	ClearSafe Comfort [®] Safety IV, 14G x 1 1/4"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84216	ClearSafe Comfort [®] Safety IV, 16G x 1 1/4"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84218	ClearSafe Comfort [®] Safety IV, 18G x 1 1/4"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84220	ClearSafe Comfort [®] Safety IV, 20G x 1 1/4"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-842202	ClearSafe Comfort [®] Safety IV, 20G x 1"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84222	ClearSafe Comfort [®] Safety IV, 22G x 1"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84224	ClearSafe Comfort [®] Safety IV, 24G x 3/4"	50 EA/BX, 4 BX/CS, 200 EA/CS