



Our I.V Administration Sets are the first in the industry to introduce a color coded packaging system, to ensure you can always find the drip rate you are looking for.




SPECIFICATION:

- Individually packaged
- Sterile
- Not made with natural rubber latex
- Not made with DEHP
- Non-pyrogenic
- Medical Grade PVC

MS-83110 - 10 Drop
MS-83115 - 15 Drop
MS-83160 - 60 Drop


50/CASE



83", Non Vented Drip Chamber, Roller Clamp, Split Septum Y-Site, Needleless PRN, Rotary Male Luer Lock

MS-831106 - 10 Drop
MS-831156 - 15 Drop
MS-831606 - 60 Drop


50/CASE



83", Non Vented Drip Chamber, Roller Clamp, Split Septum Y-Site, Needleless PRN, Rotary Male Luer Lock
With pre-attached 6" extension set

MS-83406 - 15 Drop

50/CASE




78", Non Vented Drip Chamber, Roller Clamp, Y-Site, Needleless PRN, Rotary Male Luer Lock

IV Administration Sets



MS-83401 - 10 Drop
MS-83403 - 60 Drop

50/CASE



72", Non Vented Drip Chamber, Roller Clamp, Y-Site, Needleless PRN, Rotary Male Luer Lock

MS-83171T - 60 Drop


50/CASE



88", Vented Drip Chamber, Roller Clamp, two (2) Pinch Clamps, Split Septum Y-Site, Needleless PRN, Rotary Luer Lock, Double Sided Tape

MS-83206 - 10 Drop


48/CASE



96", Non Vented Drip Chamber, Back Check Valve, Slide Clamp, Flow Regulator, two (2) Split Septum Y-Sites, Rotary Male Luer Lock

MS-83210 - 10 Drop

50/CASE



83", Non Vented Drip Chamber, Roller Clamp, two (2) Needleless PRNs, Rotary Male Luer Lock




www.medsourcelabs.com

Item #	Description	UOM
MS-83110	83", 10 Drop, Split Septum Y-Site, Needleless PRN	50 EA/CS
MS-83115	83", 15 Drop, Split Septum Y-Site, Needleless PRN	50 EA/CS
MS-83160	83", 60 Drop, Split Septum Y-Site, Needleless PRN	50 EA/CS
MS-831106	83", 10 Drop, Split Septum Y-Site, Needleless PRN, 6" ext. set	50 EA/CS
MS-831156	83", 15 Drop, Split Septum Y-Site, Needleless PRN, 6" ext. set	50 EA/CS
MS-831606	83", 60 Drop, Split Septum Y-Site, Needleless PRN, 6" ext. set	50 EA/CS
MS-83401	72", 10 Drop, Y-Site, Needleless PRN	50 EA/CS
MS-83403	72", 60 Drop, Y-Site, Needleless PRN	50 EA/CS
MS-83406	78", 15 Drop, Y-Site, Needleless PRN	50 EA/CS
MS-83171T	88", 60 Drop, Split Septum Y-Site, Needleless PRN	50 EA/CS
MS-83206	96", 10 Drop, two (2) Split Septum Y-Sites, Flow Regulator	48 EA/CS
MS-83210	83", 10 Drop, two (2) Needleless PRNs	50 EA/CS

I.V. Administration Sets, 100+ inch

MS-83175 - 10 Drop


50/CASE



100", Non Vented Drip Chamber, two (2) Slide Clamps, Roller Clamp, two (2) Split Septum Y-Sites, Male Rotary Luer Lock

MS-83105 - 15 Drop

50/CASE



105", Vented Drip Chamber, Back Check Valve, Slide Clamp, Roller Clamp, two (2) Needleless PRN Y-Sites, Male Rotary Luer Lock

MS-83122 - 15 Drop


50/CASE



122", Vented Drip Chamber, Roller Clamp, Back Check Valve, Stopcock, Needleless PRN Y-Site, Male Rotary Luer Lock

MS-83310 - 10 Drop

50/CASE



104", Vented Drip Chamber, Roller Clamp, (2) Needleless PRNs, (2) Split Septum Y-Sites, Roller Clamp, Slide Clamp, Male Rotary Luer Lock
6" Pre-Attached Extension Set

SPECIFICATION:

- Individually packaged
- Not made with natural rubber latex
- Non-pyrogenic
- Sterile
- Not made with DEHP
- Medical Grade PVC

Item #	Description	UOM
MS-83175	IV Administration Set,100"	50 EA/CS
MS-83105	IV Administration Set,105", with Back-Check Valve	50 EA/CS
MS-83122	IV Administration Set,122", with Back-Check Valve and Stopcock	50 EA/CS
MS-83310	IV Administration Set,104", with 6" Extension Set	50 EA/CS