



One-Handed Operation
Safely dispose of sharps immediately after use

Puncture-Resistant Bottom
Prevents contaminated sharps waste from contacting the user

Large Opening
Easily accommodates small syringes, IV catheter needles and other sharps waste

Small Size
Portable for EMS professionals

SPECIFICATION

- 6.7" (17cm) in Length, 1.2" (3cm) in Diameter, Empty Weight of .9oz (255grams)
- Made of Plastic
- Tested to the Standard of ISO ASTM F2132-01 for Puncture Resistance
- Locking Lid

Item #	Description	UOM
MS-64250	Sharps Dart	24/CS