


# Single Head Stethoscope



The MedSource Single Head Stethoscope is a device that does not compromise on sound quality. This unit features a polished aluminum chest piece, 22" Y-Tubing that is not made with natural rubber latex, and chrome-plated brass binaurals.



[www.medsourcelabs.com](http://www.medsourcelabs.com)

| Item #   | Description                     | UOM         |
|----------|---------------------------------|-------------|
| MS-70021 | Single Head Stethoscope, Black  | EA, 50/Case |
| MS-70022 | Single Head Stethoscope, Blue   | EA, 50/Case |
| MS-70023 | Single Head Stethoscope, Red    | EA, 50/Case |
| MS-70025 | Single Head Stethoscope, Yellow | EA, 50/Case |