

Non-Adherent
Allows the wound to
drain while letting air in



Highly Absorbant
Quickly absorbs fluids
without reopening wound

SPECIFICATION

- Various Sizes
- Highly Absorbant
- Sterilized with Ethylene Oxide
- Non-Adherent

Item #	Description	UOM
MS-GZ2212S MS-GZ4412S	Gauze, 2x2, 12 Ply Gauze, 4x4, 12 Py	50 PK/BX, 50 BX/CS, 1500 EA/CS 25 PK/BX, 50 BX/CS, 600 EA/CS