



**Blood Pressure Kit 5\***

**Carrying Case**  
Allows for easy transport

**Various Cuff Sizes**  
Fits a wide variety of patients

**Palm Aneroid Gauge**  
Easy calibration



**Blood Pressure Kit 3\***

## SPECIFICATION

- Blood Pressure Kit 5 comes with 5 cuff sizes: infant, child, adult, large adult and thigh\*
- Blood Pressure Kit 3 comes with 3 cuff sizes: child, adult and large adult\*
- Carrying case made of heavy-duty nylon
- Chrome-plated palm aneroid gauge
- One-year limited warranty

*\*Both blood pressure kits come in a carrying case with pockets for five cuffs, but Blood Pressure Kit 3 comes with three cuffs and Blood Pressure Kit 5 comes with five cuffs.*

Item #	Description	UOM
MS-MED3N	Blood Pressure Kit, 3 Cuffs, Navy	EA, 25 EA/CS
MS-MED3O	Blood Pressure Kit, 3 Cuffs, Orange	EA, 25 EA/CS
MS-MED5O	Blood Pressure Kit, 5 Cuffs, Orange	EA, 25 EA/CS
MS-MED5N	Blood Pressure Kit, 5 Cuffs, Navy	EA, 25 EA/CS