

# IV Stand User Manual

MS-IVPOLE, MS-47100 & MS-47110

MedSource 8600 Shelby Court Chanhassen, MN 55317 Made in China REV. 0

Toll Free: 1-800-876-8264 . Phone: 952-241-0131 . Fax: 952-241-0136

### 1. Applied for:


Used in hospital and clinic center.

#### 2. Technical Parameters:

- 1) With power-coated steel, stainless pole (MS-IVPOLE)
- 2) Base with ABS plastic, aluminum alloy pole (MS-47110)
- 3) Base with aluminum alloy, aluminum alloy pole (MS-47100)


# 3. Using Method:

1) The method of MS-IVPOLE



## 3. Using Method:

2) The method of MS-47100 and MS-47110





#### **Cautions:**

- 1) If damaged, repair or replace immediately.
- 2) Keep clean as routine (include disinfect).
- 3) Examine and repair the loose parts.
- 4) Store this IV Stand in a place of damp proof and non-corrosive environment.
- 5) Transfer by common transporting tools. Packed in cartons.