

Cot Sheets

Special Use



MS-44100 Poly-Tissue Disposable Sheet

SPECIFICATION

- Fluid resistant* polyethylene outer surface
- Tissue on the inside for comfort and softness
- Lightweight yet durable
- 40" x 84"



MS-42205 Heavy Duty, Cot Sheet with Pouch

SPECIFICATION

- Fluid resistant* polyethylene top layer
- Air-laid bottom layer for superior strength and support for the patient
- Elastic ends for secure fit
- Convenient pouch sewn into sheet
- 30" x 90"



MS-42210 MS-42212 Ready Fit, Cot Sheet

SPECIFICATION

- Ready fit cot sheet made especially for XPS mattress
- Non-woven polypropylene
- 36" x 90"
- MS-42210, 90 GSM
- MS-42212, 70 GSM

**Fluid could pass through the material with pressure or time*



www.medsourcelabs.com

Item #	Description	UOM
MS-44100	Poly-Tissue Flat Sheet	50 EA/CS
MS-42205	Heavy-Duty Cot Sheet with Pouch	36 EA/CS
MS-42210	Ready Fit Cot Sheet	30 EA/CS
MS-42212	Ready Fit Cot Sheet	30 EA/CS