


MS-83090 - 6"
100/CASE




Needleless PRN, Split Septum Y-Site, Slide Clamp, Male Luer Lock.
Approximate Priming Volume 1.5 mL

MS-83091 - 6"
100/CASE


MS-83092 - 8"
100/CASE

MS-83094 - 10"
100/CASE




Needleless PRN, Pinch Clamp, Rotary Male Luer Lock.
Approximate Priming Volume 2.0 mL

MS-83097 - 19"
48/CASE



Flow Regulator, Split Septum Y-Site, Slide Clamp, Rotary Male Luer Lock
Approximate Priming Volume 3.2 mL
The flow regulator dial reads from 5ml/h to 250ml/h

MS-83216 - 21"
50/CASE



Winged Male Luer, Rotary Male Luer Lock
Approximate Priming Volume 2.3 mL

SPECIFICATION

- Medical grade PVC
- Needleless

- Not made with natural rubber latex
- Not made with DEHP

- Non-pyrogenic
- Single use only

Item #	Description	UOM
MS-83090	Needleless PRN, 6", Split Septum Y Site, Side Clamp, Male Luer Lock	100/CS
MS-83091	Needleless PRN, 6", Pinch Clamp, Rotary Male Luer Lock	100/CS
MS-83092	Needleless PRN, 8", Pinch Clamp, Rotary Male Luer Lock	100/CS
MS-83094	Needleless PRN, 10", Pinch Clamp, Rotary Male Luer Lock	100/CS
MS-83097	Flow Regulator, 19", Split Septum Y Site, Side Clamp, Rotary Male Luer Lock	48/CS
MS-83216	Winged Male Luer, 21", Male Luer Lock	50/CS