


PurSafe

Safety IV Cannula with Catheter and Without Injection Valve and Without Small wings*



SPECIFICATIONS



*Not for sale in the U.S.

- PUR or PTFE catheter material
- Disposable
- Passive style safety mechanism
- Ground stainless-steel needle
- Not made with natural rubber latex
- Not made with DEHP
- Non-pyrogenic
- Sterilized with ethylene oxide
- ISO compliant color coding
- Needlestick Safety and Prevention Act (NSPA) of 2000 compliant
- Radiopaque

All catheters packaged 50/Box and 200 EA/CS


www.medsourcelabs.com

PUR MATERIAL			
Item #	Size	Item #	Size
MS-89114	14G x 1.25"	MS-891204	20G x 2.5"
MS-891142	14G x 1.77"	MS-891203	20G x 1.88"
MS-89116	16G x 1.25"	MS-891205	20G x 3"
MS-891162	16G x 1.77"	MS-89122	22G x 1"
MS-89118	18G x 1.77"	MS-891224	22G x 2.5"
MS-891182	18G x 1.25"	MS-891223	22G x 1.88"
MS-891184	18G x 2.5"	MS-891225	22G x 3"
MS-891183	18G x 3"	MS-89124	24G x 0.75"
MS-89120	20G x 1.25"	MS-891242	24G x 1.88"
MS-891202	20G x 1"	MS-891244	24G x 2.5"

PTFE MATERIAL			
Item #	Size	Item #	Size
MS-89314	14G x 1.25"	MS-893204	20G x 2.5"
MS-89314	14G x 1.77"	MS-893203	20G x 1.88"
MS-89316	16G x 1.25"	MS-893205	20G x 3"
MS-893162	16G x 1.77"	MS-89322	22G x 1"
MS-89318	18G x 1.77"	MS-893224	22G x 2.5"
MS-893182	18G x 1.25"	MS-893223	22G x 1.88"
MS-893184	18G x 2.5"	MS-893225	22G x 3"
MS-893183	18G x 3"	MS-89324	24G x 0.75"
MS-89320	20G x 1.25"	MS-893242	24G x 1.88"
MS-893202	20G x 1"	MS-893244	24G x 2.5"

PurSafe

Safety IV Cannula with Catheter and With Injection Port and Small wings*



Passive Fail-Safe Safety Clip
Automatically activated as needle is removed from catheter

Flashback Chamber
Allows for clear visual confirmation of successful venipuncture

Injection Port
Color-coded cap port

Removable Flashplug

Triple Facet, Back Cut Beveled Needle
Designed to provide ease of insertion

Thumb Stopper
Helps with one-handed insertion techniques

Gently Tapered Catheter Tip
Provides smooth transition from needle to catheter during insertion

Suture Wings
Helps to keep catheter stable**



SPECIFICATIONS

*Not for sale in the U.S.
**When used with securement device

- PUR or PTFE catheter material
- Disposable
- Sterilized with ethylene oxide
- Not made with natural rubber latex
- Not made with DEHP
- Non-pyrogenic
- Ground stainless-steel needle
- Passive style safety mechanism
- Radiopaque
- Winged hub with injection port
- ISO compliant color coding
- Needlestick Safety and Prevention Act (NSPA) of 2000 compliant

All catheters packaged 50/Box and 200 EA/CS

www.medsourcelabs.com

PUR MATERIAL			
Item #	Size	Item #	Size
MS-89214	14G x 1.25"	MS-892204	20G x 2.5"
MS-89214	14G x 1.77"	MS-892203	20G x 1.88"
MS-89216	16G x 1.25"	MS-892205	20G x 3"
MS-892162	16G x 1.77"	MS-89222	22G x 1"
MS-89218	18G x 1.77"	MS-892224	22G x 2.5"
MS-892182	18G x 1.25"	MS-892223	22G x 1.88"
MS-892184	18G x 2.5"	MS-892225	22G x 3"
MS-892184	18G x 3"	MS-89224	24G x 0.75"
MS-89220	20G x 1.25"	MS-892242	24G x 1.88"
MS-892202	20G x 1"	MS-892244	24G x 2.5"

PTFE MATERIAL			
Item #	Size	Item #	Size
MS-89414	14G x 1.25"	MS-892204	20G x 2.5"
MS-894142	14G x 1.77"	MS-892203	20G x 1.88"
MS-89416	16G x 1.25"	MS-892205	20G x 3"
MS-894162	16G x 1.77"	MS-89222	22G x 1"
MS-89418	18G x 1.77"	MS-892224	22G x 2.5"
MS-894182	18G x 1.25"	MS-892223	22G x 1.88"
MS-892184	18G x 2.5"	MS-892225	22G x 3"
MS-892183	18G x 3"	MS-89224	24G x 0.75"
MS-89220	20G x 1.25"	MS-892242	24G x 1.88"
MS-892202	20G x 1"	MS-892244	24G x 2.5"