


ClearSafe Comfort® BC



Slide Style Safety IV Catheter

Blood Control with Active Check Valve Technology



SPECIFICATIONS



- Incorporated blood control
- PUR catheter material
- Slide style needle retraction
- Ground stainless-steel needle
- Straight hub
- Disposable
- Non-pyrogenic
- Not made with natural rubber latex
- Not made with DEHP
- ISO compliant color coding
- Needlestick Safety and Prevention Act (NSPA) of 2000 compliant



www.medsourcelabs.com

Item #	Description	UOM
MS-86214	ClearSafe Comfort® Blood Control Safety IV, 14G x 1 1/4"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-86216	ClearSafe Comfort® Blood Control Safety IV, 16G x 1 1/4"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-86218	ClearSafe Comfort® Blood Control Safety IV, 18G x 1 1/4"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-86220	ClearSafe Comfort® Blood Control Safety IV, 20G x 1 1/4"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-862202	ClearSafe Comfort® Blood Control Safety IV, 20G x 1"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-86222	ClearSafe Comfort® Blood Control Safety IV, 22G x 1"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-86224	ClearSafe Comfort® Blood Control Safety IV, 24G x 3/4"	50 EA/BX, 4 BX/CS, 200 EA/CS