


Over-the-Ear Design
Maintains proper positioning
of nasal tips

Star Lumen Supply Tubing
Durable and kink-resistant



SPECIFICATION

- Various Sizes
- For Patients Requiring Medium or Low Concentrations of Oxygen
- Oxygen Concentration Delivery: 24-44%
- Recommended Flow Rate: 1-6 LPM

Item #	Description	UOM
MS-24003	Nasal Cannula, Non-Flared Tip, 7'	50 EA/CS
MS-24004	Nasal Cannula, Flared Tip, 7'	50 EA/CS
MS-24005	Nasal Cannula, Non-Flared Tip, 14'	50 EA/CS
MS-24008	Nasal Cannula, Non-Flared Tip w/ Universal Connector, 7'	50 EA/CS
MS-24101	Nasal Cannula, Pediatric, 7'	50 EA/CS
MS-20007	Star Lumen O2 Tubing Only, 7'	50 EA/CS