



Large LCD Screen
Provides clear reading of results



Glucose Test Strips
Displays blood test results in six seconds

SPECIFICATION

- Dimensions: 94.5mm x 56mm x 27.5mm
- Weight: 56g
- Display: large LCD
- Auto electrode inserting detection
- Auto turn-off after one minute without action
- Comes with one CR2023 battery
- Glucometer comes in a hard case that measures 6" L x 3" W x 1.5" H

*Items sold separately.
Please consult instructions for use.*



www.medsourcelabs.com

Item #	Description	UOM
MS-76010	Glucometer	EA
MS-76012	Glucometer Test Strips	25 Strips/Vials, 2 Vials/BX
MS-76014	Glucometer Lancet	100 EA/BX
MS-76016	Glucometer Test Solution	EA