



Rip-Stop Fabric
Ensures longevity

Color Coded Straps
Ease of application

Also comes in red.

SPECIFICATION

- Color coded straps
- Can be used on pregnant and pediatric patients
- Includes head support, reusable head straps, and a carrying bag
- Internal lining: high tenacity polyester, 205 GSM
- Outer cover: PVC two-sided coating, 380 GSM, polyester, waterproof at 220 PSI: no leak
- Lifting handles: 2" polyester webbing belt
- Spine cushion: PVC two-sided coating with U foam insert (30 density)
- Wood work: 30 sticks
 - Main frame: teak wood
 - Outer support: pine wood
- Carton weight: 15.2 KG

Item #	Description	UOM
MS-ED2253	Extrication Device, Green	4 EA/CS
MS-ED3000	Extrication Device, Red	4 EA/CS