



SPECIFICATION

- High Strength Aluminum Alloy 6061-T6
- Standard Finish: Brushed Body with a Clear Coat Finish
- Fully DOT Approved

	CAPACITY/Cu. ft/Liters	Length	Diameter	Empty Weight Lb/Kg	Thread Size
MS-MED-ET:	24.0/679	25.4/645	4.38/111.3	7.8/3.6	0.750 - 16UNF - 2B
MS-MED-DT:	15.0/425	16.7/424	4.38/111.3	5.4/2.4	0.750 - 16UNF - 2B

Item #	Description	UOM
MS-MED-DT	Cylinder, D-Tank w/ Toggle	EA
MS-MED-ET	Cylinder, E-Tank w/ Toggle	EA