


Safety Tip
Prevents injury to the patient

One Serrated Edge
Cuts through a variety of material



Molded Plastic Handle
Provides comfort for the user

SPECIFICATION

- 7 1/4" in length
- Hardened stainless-steel blades
- Available in nine assorted colors
- Autoclavable

Item #	Description	UOM
MS-SH001B	EMT Shears, Black	200/CS
MS-SH002BL	EMT Shears, Blue	200/CS
MS-SH003P	EMT Shears, Hot Pink	200/CS
MS-SH004Y	EMT Shears, Yellow	200/CS
MS-SH005R	EMT Shears, Red	200/CS
MS-SH006G	EMT Shears, Green	200/CS
MS-SH007PR	EMT Shears, Purple	200/CS
MS-SH008O	EMT Shears, Orange	200/CS
MS-SH009W	EMT Shears, White	200/CS