


Uncuffed and Cuffed with Stylet

MedSource Endotracheal Tubes are engineered with clear, flexible PVC, and a smooth Murphy eye as an alternative vent port



SPECIFICATION

- Standard 15mm connector
- Oral or nasal use
- Short or long-term use
- Medical grade PVC
- Sterilized for single-use
- Not made with natural rubber latex



www.medsourcelabs.com

Item No.	ET Tube Uncuffed w/Stylet Size	UOM	Item No.	ET Tube Cuffed w/Stylet Size
MS-23325	2.5mm	10 EA/BX	MS-23450	5.0mm
MS-23330	3.0mm	10 BX/CS	MS-23455	5.5mm
MS-23335	3.5mm	100 EA/CS	MS-23460	6.0mm
MS-23340	4.0mm		MS-23465	6.5mm
MS-23345	4.5mm		MS-23470	7.0mm
MS-23350	5.0mm		MS-23475	7.5mm
MS-23355	5.5mm		MS-23480	8.0mm
			MS-23485	8.5mm
			MS-23490	9.0mm
			MS-23495	9.5mm
			MS-23410	10.0mm