

Single-Use Emergency Blankets



Reflective Surface
Provides high visibility

Mylar



Poly-Foam
Designed to keep the patient warm and dry

Poly-Foam

Water Resistant* Backing
Protects the patient from liquids



Poly-Tissue

SPECIFICATION

- MS-B100: 52" x 83"
- MS-B200: 58" x 90"
- MS-B300: 56" x 90"
- Three Types: Mylar, Poly-Foam and Poly-Tissue

**Fluid may pass through material with pressure or time*



www.medsourcelabs.com

Item #	Description	UOM
MS-B100	Mylar Emergency Blanket	200 EA/CS
MS-B200	Poly-Foam Emergency Blanket, Yellow	18 EA/CS
MS-B300	Poly-Tissue Emergency Blanket, Yellow	24 EA/CS