

One Quart Sharps Container



Transparent Locking Lid

Lid locks to secure biohazard waste of needles, syringes and IV catheters

Translucent Red ABS Plastic With Fill level Identification

Puncture and leak resistant material is also transparent for fill level identification

One Quart Size

Compact size fits easily into gear bags, but large enough to use multiple times over

SPECIFICATION

Materials	ABS Plastic
Dimensions	Height: 6.25" Depth: 4.25" Width: 4.25"
Product Testing	Tested to the Standard of ISO ASTM F2132-01 for Puncture Resistance
ISO Standards	Meets ISO 23907:2012 Sharps Protection Standard



www.medsourcelabs.com

Item #	Description	UOM
MS-64251	One Quart Sharps Container	60 EA/CS