


ClearSafe Comfort[®] BC



Slide Style Safety IV Catheter

Blood Control with Active Check Valve Technology



SPECIFICATIONS



- Incorporated blood control
- PUR catheter material
- Slide style needle retraction
- Ground stainless-steel needle
- Straight hub
- Disposable
- Non-pyrogenic
- Not made with natural rubber latex
- Not made with DEHP
- ISO compliant color coding
- Needlestick Safety and Prevention Act (NSPA) of 2000 compliant
- Radiopaque



www.medsourcelabs.com

Item #	Gauge and Length	Catheter ID (mm)	Catheter OD (mm)	Gravity Flow Rate (mL/min)	UOM
MS-86214	14G x 1.25"	1.6 mm	2.1 mm	297 mL/min	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-86216	16G x 1.25"	1.38 mm	1.8 mm	220 mL/min	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-86218	18G x 1.25"	0.95 mm	1.3 mm	94 mL/min	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-86220	20G x 1.25"	0.78 mm	1.1 mm	61 mL/min	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-862202	20G x 1"	0.78 mm	1.1mm	61 mL/min	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-86222	22G x 1"	0.62 mm	0.9 mm	36 mL/min	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS
MS-86224	24G x 0.75"	0.49 mm	0.7 mm	20 mL/min	EA, 50 EA/BX, 4 BX/CS, 200 EA/CS