


ClearSafe® Safety IV Catheter



SPECIFICATION


- PTFE catheter material
- Slide style needle retraction
- Radiopaque
- Non-pyrogenic, not made with DEHP
- Not made with natural rubber latex
- Disposable
- Siliconized ground stainless steel needle
- Straight hub
- ISO compliant color coding
- Needle Stick Prevention and Safety Act (NSPA) of 2000 compliant



www.medsourcelabs.com

| Item # | Needle Gauge and Length | Catheter ID (mm) | Catheter OD (mm) | Gravity flow rate (mL/min) | U of M |
|-----------|-------------------------|------------------|------------------|----------------------------|----------------------------|
| MS-84114 | 14G x 1.25" | 32 mm | 2.1 mm | 297 mL/min | 50/BX, 200/CS, (4 BX/Case) |
| MS-84116 | 16G x 1.25" | 32 mm | 1.8 mm | 220 mL/min | 50/BX, 200/CS, (4 BX/Case) |
| MS-84118 | 18G x 1.25" | 32 mm | 1.3 mm | 94 mL/min | 50/BX, 200/CS, (4 BX/Case) |
| MS-84120 | 20G x 1.25" | 32 mm | 1.1 mm | 61 mL/min | 50/BX, 200/CS, (4 BX/Case) |
| MS-841202 | 20G x 1" | 25 mm | 1.1 mm | 61 mL/min | 50/BX, 200/CS, (4 BX/Case) |
| MS-84122 | 22G x 1" | 25 mm | 0.9 mm | 36 mL/min | 50/BX, 200/CS, (4 BX/Case) |
| MS-84124 | 24G x 0.75" | 19 mm | 0.7 mm | 20 mL/min | 50/BX, 200/CS, (4 BX/Case) |

ClearSafe Comfort® Safety IV Catheter



SPECIFICATION

- PUR catheter material
- Slide style needle retraction
- Radiopaque
- Non-pyrogenic, not made with DEHP
- Not made with natural rubber latex
- Disposable
- Siliconized ground stainless steel needle
- Straight hub
- ISO compliant color coding
- Needle Stick Prevention and Safety Act (NSPA) of 2000 compliant



www.medsourcelabs.com

| Item # | Needle Gauge and Length | Catheter ID (mm) | Catheter OD (mm) | Gravity flow rate (mL/min) | U of M |
|-----------|-------------------------|------------------|------------------|----------------------------|----------------------------|
| MS-84214 | 14G x 1.25" | 32 mm | 2.1 mm | 297 mL/min | 50/BX, 200/CS, (4 BX/Case) |
| MS-84216 | 16G x 1.25" | 32 mm | 1.8 mm | 220 mL/min | 50/BX, 200/CS, (4 BX/Case) |
| MS-84218 | 18G x 1.25" | 32 mm | 1.3 mm | 94 mL/min | 50/BX, 200/CS, (4 BX/Case) |
| MS-84220 | 20G x 1.25" | 32 mm | 1.1 mm | 61 mL/min | 50/BX, 200/CS, (4 BX/Case) |
| MS-842202 | 20G x 1" | 25 mm | 1.1 mm | 61 mL/min | 50/BX, 200/CS, (4 BX/Case) |
| MS-84222 | 22G x 1" | 25 mm | 0.9 mm | 36 mL/min | 50/BX, 200/CS, (4 BX/Case) |
| MS-84224 | 24G x 0.75" | 19 mm | 0.7 mm | 20 mL/min | 50/BX, 200/CS, (4 BX/Case) |