


ClearSafe Comfort® CS

Slide Style Closed System
Safety IV Catheter with Dual Port



SPECIFICATION



**when used with securement device*

- | | | |
|---|--|--|
| <ul style="list-style-type: none">• Pre-attached 6 in. IV extension set• Needleless PRN• PUR catheter material• Slide style needle retraction• Radiopaque | <ul style="list-style-type: none">• Non-pyrogenic, not made with DEHP• Not made with natural rubber latex• Disposable• Siliconized ground stainless steel needle• Winged hub | <ul style="list-style-type: none">• ISO compliant color coding• Needle Stick Prevention and Safety Act (NSPA) of 2000 compliant |
|---|--|--|



www.medsourcelabs.com

Item #	Gauge size	Length in mm	Length in inches	Catheter ID (mm)	Catheter OD (mm)	Gravity Flow Rate (mL/min)	UOM
MS-87318	18G	32 mm	1.25"	0.95 mm	1.3 mm	94 mL/min	EA, 20 EA/BX, 4 BX/CS, 80 EA/CS
MS-87320	20G	32 mm	1.25"	0.78 mm	1.1 mm	61 mL/min	EA, 20 EA/BX, 4 BX/CS, 80 EA/CS
MS-87322	22G	25 mm	1"	0.62 mm	0.9 mm	36 mL/min	EA, 20 EA/BX, 4 BX/CS, 80 EA/CS
MS-87324	24G	19 mm	0.75"	0.49 mm	0.7 mm	20 mL/min	EA, 20 EA/BX, 4 BX/CS, 80 EA/CS