


ClearSafe Comfort® CS



Slide Style Closed System Safety IV Catheter with Single Port



SPECIFICATION



*when used with a securement device

- Pre-attached 6 in. IV extension set
 - Needle-free connector
 - PUR catheter material
 - Slide style needle retraction
 - Radiopaque
- Non-pyrogenic, not made with DEHP
 - Not made with natural rubber latex
 - Disposable
 - Siliconized ground stainless steel needle
 - Winged hub
- ISO compliant color coding
 - Needle Stick Prevention and Safety Act (NSPA) of 2000 compliant




www.medsourcelabs.com

Item #	Gauge size	Length in mm	Length in inches	Catheter ID (mm)	Catheter OD (mm)	Gravity Flow Rate (mL/min)	UOM
MS-87218	18G	32 mm	1.25"	0.95 mm	1.3 mm	94 mL/min	EA, 20 EA/BX, 4 BX/CS, 80 EA/CS
MS-872202	20G	25 mm	1"	0.78 mm	1.1 mm	61 mL/min	EA, 20 EA/BX, 4 BX/CS, 80 EA/CS
MS-87220	20G	32 mm	1.25"	0.78 mm	1.1 mm	61 mL/min	EA, 20 EA/BX, 4 BX/CS, 80 EA/CS
MS-87222	22G	25 mm	1"	0.62 mm	0.9 mm	36 mL/min	EA, 20 EA/BX, 4 BX/CS, 80 EA/CS
MS-87224	24G	19 mm	0.75"	0.49 mm	0.7 mm	20 mL/min	EA, 20 EA/BX, 4 BX/CS, 80 EA/CS

ClearSafe Comfort® CS

Slide Style Closed System
Safety IV Catheter with Dual Port

MEDSOURCE



SPECIFICATION



*when used with securement device

- Pre-attached 6 in. IV extension set
- Needleless PRN
- PUR catheter material
- Slide style needle retraction
- Radiopaque
- Non-pyrogenic, not made with DEHP
- Not made with natural rubber latex
- Disposable
- Siliconized ground stainless steel needle
- Winged hub
- ISO compliant color coding
- Needle Stick Prevention and Safety Act (NSPA) of 2000 compliant

MEDSOURCELABS

www.medsourcelabs.com

Item #	Gauge size	Length in mm	Length in inches	Catheter ID (mm)	Catheter OD (mm)	Gravity Flow Rate (mL/min)	UOM
MS-87318	18G	32 mm	1.25"	0.95 mm	1.3 mm	94 mL/min	EA, 20 EA/BX, 4 BX/CS, 80 EA/CS
MS-87320	20G	32 mm	1.25"	0.78 mm	1.1 mm	61 mL/min	EA, 20 EA/BX, 4 BX/CS, 80 EA/CS
MS-87322	22G	25 mm	1"	0.62 mm	0.9 mm	36 mL/min	EA, 20 EA/BX, 4 BX/CS, 80 EA/CS
MS-87324	24G	19 mm	0.75"	0.49 mm	0.7 mm	20 mL/min	EA, 20 EA/BX, 4 BX/CS, 80 EA/CS