



The MedSource CPR Board is designed for CPR training and administration. It is anatomically contoured to help lift and position the patient for effective CPR.



### Instructions:

- Gently raise patients head.
- Slide CPR Board under patients head.
- Commence CPR on patient.

Manual for Item Number(s):

MS-CPRB CPR Board

MedSource Europe:  
Robert-Bosch-str 19, 73117  
Wangen, Germany  
Tel: +49 71 61 - 98 34 999

MedSource USA:  
8600 Shelby Court  
Chanhassen MN, 55317 USA  
Tel: 1-952-472-0131  
export@medsourcelabs.com  
Made in China

MedSource UAE:  
Le Solarium Centre, Office 916,  
Dubai Silicon Oasis, Dubai, UAE  
Tel: 971-4-320 7123