



IV Stand User Manual

MS-IVPOLE, MS-47100 & MS-47110



MedSource Labs, LLC
8600 Shelby Court
Chanhassen MN, 55317 USA
MSL-REG-LQR-10082, REV. 02
Toll Free: 1-800-876-8264
Phone: 952-241-0131
Fax: 952-241-0136



Made in China

1. Applied for:


Used in hospital and clinic center.

2. Technical Parameters:

- 1) With power-coated steel, stainless pole (MS-IVPOLE)
- 2) Base with ABS plastic, aluminum alloy pole (MS-47110)
- 3) Base with aluminum alloy, aluminum alloy pole (MS-47100)


3. Using Method:

- 1) The method of MS-IVPOLE



3. Using Method:

2) The method of MS-47100 and MS-47110



Cautions:

- 1) If damaged, repair or replace immediately.
 - 2) Keep clean as routine (include disinfect).
 - 3) Examine and repair the loose parts.
 - 4) Store this IV Stand in a place of damp proof and non-corrosive environment.
 - 5) Transfer by common transporting tools. Packed in cartons.
- Shelf Life: 8 years