

Cot Sheets

Heavy Duty



MS-42604 Pro-fitted Cot Sheet, Heavy Duty, Sewn Ends

SPECIFICATION

- Made from a strong and durable emergency linen, 45 gsm
- Provides fluid resistant* barrier
- Sewn corners for a secure fit



MS-001HS Cot Sheet, Heavy Duty, Square Ends

SPECIFICATION

- Made from a spun bonded polypropylene, 95 gsm
- Fluid resistant*
- Specifically designed for emergency service stretchers

**Fluid could pass through the material with pressure or time*



www.medsourcelabs.com

Item #	Description	UOM
MS-42604	Cot Sheet, Heavy Duty, Fitted, Sewn Ends, 30" x 72"	50 EA/CS
MS-001HS	Cot Sheet, Heavy Duty, Fitted, Square Ends, 30" x 72"	50 EA/CS