


NIK Kit

Needle Decompression Indicator Kit



Buy American Act Compliant- Manufactured in accordance with federal sourcing requirements for domestic procurement

SPECIFICATIONS

- Needle Size: 14, 12, 10 Gauge x 3.25 inch
- PTFE Catheter material
- EO Sterilization
- Single Use Only
- Not made with natural rubber latex
- Not made with DEHP



www.medsourcelabs.com

Item #	Description	UOM
MS-87110	NIK Kit - 10G	EA
MS-87112	NIK Kit - 12G	EA
MS-87114	NIK Kit - 14G	EA