




SPECIFICATIONS

- Lightweight One-Piece Molded Board
- Tough, Space-Age Materials
- Floats for Water Rescue
- 100% X-Ray Translucent
- Load limit: 600 lbs
- Length of Spineboard: 72"

Item #	Description	UOM
MS-93000	Spineboard, Green, 16" W	EA
MS-93000-O	Spineboard, Orange, 16" W	EA
MS-93000Y-16	Spineboard, Yellow, 16" W	EA
MS-93000-18O	Spineboard, Orange, 18" W	EA
MS-93000-18Y	Spineboard, Yellow, 18" W	EA